

Distributed by: Medset Specialties Ltd.

**NEXT
GENERATION**

FlowTrax® FT-2 IV Pump Analyzer

NEW

**Electronic Occlusion, Quick Prime
& Field-Serviceable Maintenance Features**



Infusion Testing

Results within 1% accuracy in 3 minutes or less

NEW Electronic Occlusion feature (0-40 PSI)

*NEW Faster priming with no stopping
to re-prime between tests!*

Full-Featured Pressure Meter

*Standalone pressure meter: test any type
of medical device for pressure or vacuum*

Test fluid, gas or air pressure (-12.5 to +75 PSI)

YSI 400 Temperature Meter

Measure fluid, gas or air (-28.9 to +115.5 °C)

Digital Timer and Stop Watch

*Automatically starts & stops
with fluid measurements*

All features are NIST Traceable.

**Multifunction IV Analyzer —
More Versatility, Less Cost**

PRONK TECHNOLOGIES

<http://medset.ca>

sales@medset.ca

MANUFACTURED IN U.S.A.