

NEW

SimSlim™ Multi-Parameter Simulator



10 Year Battery Life!

Revolutionizing biomedical test equipment

Introducing the SimSlim Multi-parameter Simulator:
One Button. Fifty Simulations.


- ▶ 12 Lead **ECG/Resp** with built in banana plug interface
- ▶ 4 Channels of Invasive **Blood Pressure**
- ▶ **Cardiac Output** with no adapters needed
- ▶ YSI 400/700 **Temperature**
- ▶ **Pocket sized**, rugged, reliable and **easy to use**
- ▶ **10 Year Battery** Life* (no charging or changing)

*10 years is approximate battery life depending on usage.

One button, 50 possibilities

The **SimSlim** is a compact, rugged, easy to use multi-parameter patient simulator combining all of the following features:

- **12 Lead ECG:** 8 selectable Heart Rates, 8 selectable arrhythmias
- **Respiration:** 8 selectable respiration rates
- **4 Invasive Blood Pressure Channels:** Dynamic: 120/80, 60/30, 30/10, 20/10, 125/15, 60/5, 8/2 and 4/0. Static pressures: 0, 25, 50, 100, 150, 200, 250
- **Cardiac Output Simulation:** 3, 5 and 7 (l/min)
- **Temp Simulations- YSI 400 and 700:** 25, 37 and 40 degrees C
- **10 Year Battery Life** (no charging or changing required)
- Simple to use, single button operation
- 4mm and 1/8in banana plug grooves for easy hook up to 12 lead EKG carts
- 50 different selectable simulations



OxSim™
SpO2 Simulator

SimSlim™
Multi-Parameter Simulator



SimCube™
NIBP Simulator

Combine **SimSlim** with **SimCube** and **OxSim**, save money and take your test bench with you wherever you go.



SC5+OxSlimKit



NIBP • 12 Lead ECG • Resp • 4 IBP Channels • SpO2 • Arrhythmia • Cardiac Output • Temp • Manometer



SC5+OxSlimKit Includes:
SC-5, OxSim, SimSlim, Battery Boost, carrying case and all NIBP Adapters.



Leave the Cart Behind...

PRONK TECHNOLOGIES
Medset Specialties Ltd.

Tel: 1-800-331-3187

3334 Mainway Drive, Burlington, ON L7M 1A7

Fax: 905-335-3356