

Biomedical Batteries

Medset Specialties Ltd.
3334 Mainway Drive
Burlington, ON L7M 1A7

Tel: 800-331-3187
Fax: 905-335-4456

Health Canada Establishment License # 1810

Email: sales@medset.ca
Web: <http://medset.ca>

COMPANY	MODEL	VOLTS	AH	BATTERY
Accutome	00-1912	0.00	0.00	OEM ONLY
	24-5100	3.60		LS14250BA
	PachPen	3.60		LS14250BA
Accutor	9500	9.60	2.00	MS8200
Acme Medical	All Scales *	12.00	1.20	L12V1.2F1-04
Aequitron	9550	9.60	2.30	MSL-9550
	Suction 8850	12.00	3.40	L12V3.4F1-04
	Vents LP3, 4	12.00	2.50	GL8108
Aerogen	Aeroneb	4.80	0.45	MSL-AERONEB
Aeros	Aspirator	12.00	3.40	L12V3.4F1-04
	Tote-L-Vac (4900)	12.00	2.00	L12V2.0C2-17
AeroScout	Asset RFID Tag	3.60		LS14250BA
	T2S Senso Tag	0.00	0.00	OEM ONLY
Air Shields	101	4.80	0.20	MS101
	C100 AlarmBU	7.20	0.15	M9
	C450	0.00	0.00	OEM ONLY
	C-450	0.00	0.00	OEM ONLY
	JM102	0.00	0.00	OEM ONLY
	System 5	6.00	3.40	L06V3.4F1-01
	T1-100 Incub	12.00	26.00	L12V26B1-16
	TI 500	12.00	26.00	L12V26B1-16
Airborne	20H, Voyager	12.00	26.00	MSL-5970
	Isolette 185, 185A	12.00	28.00	L12V28B1-15
	Isolette 185A Plus	12.00	26.00	MSL-5970
	T201	12.00	13.00	MSL-5722
Aitecs	2015	0.00	0.00	OEM ONLY
Alaris/Ivac	2000	4.80	0.50	MSL-AA/PX4
	4200	12.00	2.50	MSL-5167
	4400	6.00	5.00	L06V4.5F1-01
	4410	6.00	5.00	L06V4.5F1-01

* May require multiple batteries

COMPANY	MODEL	VOLTS	AH	BATTERY
Alaris/Ivac				
	4415	6.00	5.00	L06V4.5F1-01
	530	12.00	0.50	M12/500
	560	12.00	2.50	GL8108-P
	710	8.00	2.50	GL825P
	7100/7230	12.00	1.80	MSL-5860
	7200	12.00	1.80	MSL-5860
	8000 MEDLEY	12.00	4.50	MSL-5073
	821	4.80	0.60	MB5076
	Express 3	3.00	0.35	LTC-3PN
	P2000	6.00	3.20	L06V3.2F1-06
	P3000	6.00	3.00	L06V3.0F1-02
	P4000	6.00	3.20	L06V3.2F1-03
	P5000	6.00	3.00	L06V3.0F1-02
	P7000	6.00	3.40	L06V3.4F1-01
	PCAM	6.00	3.40	L06V3.4F1-01
	Signature	12.00	1.80	MSL-5860
	Syringe 700	8.00	2.50	MSL-5016
Alkaline				
	9 Volt	9.00	0.00	MSL-ALK-9V
	AA	1.50	0.00	MSL-ALK-AA
	AAA	1.50	0.00	MSL-ALK-AAA
	C	1.50	0.00	MSL-ALK-C
	D	1.50	0.00	MSL-ALK-D
	N	1.50	0.00	MSL-ALK-N
Allied Healthcare				
	L190	12.00	3.40	L12V3.4F1-04
	L190-GR	12.00	3.40	L12V3.4F1-04
American Optical				
	Pulsar IV	12.00	2.20	M0923
Amsco				
	3080/3085 *	12.00	26.00	L12V26B1-16
Arjo				
	*	12.00	4.50	L12V5.0F1-03
	*	12.00	4.00	L12V4.1F1-02
	HMX27			MSL-PL-AB-00600
	HMX446			MSL-PL-AB-00700
	KKA1100-04	24.00	4.00	MSL-PL-AB-00400
	KPA0100			MSL-PL-AB-00200
	KTA0102	24.00	2.30	MSL-PL-AB-00500
	MSP-A-AH-01010*	12.00	2.10	L12V2.3C1-17
	NDA0100-20	24.00	4.00	MSL-PL-AB-00100
	NDA0200	24.00	5.00	MSL-PL-AB-00120
	NDA0200	24.00	5.50	MSL-PL-AB-00110
	NEA0100-083			MSL-PL-AB-00300
	Opera	6.00	5.00	L06V4.5F1-01
	Opera (req 4)	6.00	5.00	0809-0012
	Sara	24.00	4.00	MSL-PL-AB-00400
	Sara Plus			MSL-PL-AB-00200

COMPANY	MODEL	VOLTS	AH	BATTERY
Armstrong	Duet	12.00	2.80	MSL-5052
	Suction Pump	12.00	4.50	L12V5.0F1-03
Array Medical	A2P	0.00	0.00	OEM ONLY
Arrow	ACADI	12.00	18.00	L12V18B1-14
	Autocat 2 Wave IABP	12.00	18.00	L12V18B1-14
Aspect	BIS Vista	0.00	0.00	OEM ONLY
	BISA2000	7.20	2.10	MSL-A2000
B. Braun	Horizon	12.00	2.00	L12V2.0F1-08
	Nxt	12.00	2.00	L12V2.0F1-08
	Nxt B/U Alarm	0.00	0.00	OEM ONLY
	Nxt B/U Mem	3.00	0.30	MSL-5739
	Outlook 100, 200	12.00	2.60	L12V2.6F1-02
	Outlook 300ES	0.00	0.00	OEM ONLY
	Outlook 400ES	0.00	0.00	OEM ONLY
	Perfusor Compact	1.20	0.00	MSL-AA
	Perfusor Compact NIC			N/A
	Space IVPump	4.80	2.10	MSL-8713180
	Vista Basic	0.00	0.00	OEM ONLY
Bard	Bladder Scan 500	6.00	5.00	0809-0012
	Bladder Scan BVI2500	6.00	2.50	0819-0012
	Harvard PCA	12.00	1.20	L12V1.2F1-04
	PCA2	4.80	5.00	MSL-PCA2
Baxter	2015 Aitecs	0.00	0.00	OEM ONLY
	6100	6.00	12.00	L06V12F1-02
	6200*	2.00	6.00	PS-260
	6201	12.00	2.00	L12V2.0F1-08
	6300	12.00	3.40	L12V3.4F1-04
	6301	12.00	3.40	L12V3.4F1-04
	6301/BU	0.00	0.00	OEM ONLY
	Aitecs 2015	0.00	0.00	OEM ONLY
	AS20-GH	6.00	0.60	MSL-5436
	AS40A	6.00	1.20	MSL-AS40A
	AS50A	6.00	1.20	MSL-AS40A
	Colleague *	12.00	2.00	L12V2.0F1-08
	Colleague B/U	0.00	0.00	OEM ONLY
	OxiSat SM020	6.00	5.00	L06V4.5F1-01
	PCA Pump	4.80	5.00	MSL-PCA2
	PCA2	4.80	5.00	MSL-PCA2
	S00225	7.20	2.00	MSL-S00225-I
	SEP-***	9.60	1.30	MSL-6307
	Sigma Spectrum	7.20	2.00	MSL-S00225-I
BCI (Biochem)				

COMPANY	MODEL	VOLTS	AH	BATTERY
BCI (Biochem)				
	3202 3302	6.00	1.30	MSL-5852
	3302	6.00	1.30	L06V1.3F1-09
	6004	6.00	3.00	MSL-5945
	8200	9.60	1.60	MSL-5903
	Arcona	0.00	0.00	MSL-AAPLQ
	Capnocheck Plus	6.00	4.40	MSL-5108
	Mini Toor Plus	6.00	3.00	MSL-5945
Bear				
	Bear II Vent	4.80	0.60	MSL-5162
	Cub 1580	4.80	0.60	MSL-5162
BHM Medical				
	550	12.00	7.20	L12V7.2F1-04
	Liberte 2	12.00	7.20	L12V7.2F1-04
	V4 *	12.00	4.50	L12V5.0F1-03
	Voyageur 4vi	12.00	4.50	L12V5.0F1-03
Biomedical Devic				
	Crossvent 3	0.00	0.00	OEM ONLY
	Crossvent 3+ (3 wire)	0.00	0.00	OEM ONLY
	Crossvent 4, 4+	0.00	0.00	OEM ONLY
Bio-Tek				
	501	9.60	1.30	MSL-5645
Bird				
	T-Bird *	12.00	0.80	MSL-5949
BpTRU Medical				
	BPM-200	8.00	2.70	L08V2.5F1-04
	BPM-300	8.00	2.70	L08V2.5F1-04
Burdick				
	ATRIA 3000	14.40	1.60	MSL-6053
	E560	0.00	0.00	OEM ONLY
	E600 E500	9.60	0.60	MSL-5565
	Eclipse 8	16.80	1.30	MSL-5876
	Eclipse 850	16.80	0.70	MSL-5834
	Eclipse Plus	16.80	2.30	MSL-5891
	Eclipse Plus Premier	16.80	4.40	MSL-5003
	Elite II	12.00	1.00	LP6/200
	Elite II / EK10	12.00	1.80	MSL-5226
	Medic 4	24.00	2.50	MSL-5215
	Medic 5	0.00	0.00	OEM ONLY
	Vision 5L	3.00		MSL-VISION5L
Cardiotronics				
		3.60	0.00	MSL-3/V100R/T
CAS Medical				
	730, 740	7.20	4.00	MSL-6036
Centron				
	BN-122	12.00	2.00	MSL-BN122
Chattanooga				
	Laser	12.00	2.60	L12V2.6F1-04
Cio-Teck				

COMPANY	MODEL	VOLTS	AH	BATTERY
Cio-Teck				
	918850	12.00	3.40	L12V3.4F1-04
Clinitec				
	Feeding Pmp	12.00	1.20	L12V1.2F1-04
Colin				
	8800 Clock	0.00	0.00	OEM ONLY
	BP8800	12.00	2.60	L12V2.6F1-02
	BP880NXT	12.00	2.00	L12V2.0C1-17
	PressMate	12.00	2.60	L12V2.6F1-02
	Prodigy	6.00	5.00	MSL-6014
Concord				
	Patient Lift (Req 2)	6.00	5.00	L06V4.5F1-01
Cook Medical				
	DP-M350 Doppler	0.00	0.00	OEM ONLY
Cordis				
	Pacer Pgrm	16.80	1.20	MSL-255A
Corometrics				
	115	4.00	0.50	L04V1/2F1-01
	145	3.50	0.70	90-M35V750
	220 pH Meter	3.50	0.70	MSL-5359
	332	7.20	0.60	MSL-5729
	502	6.00	2.50	MSL-5356
	511	6.00	7.20	L06V7.2F1-02
Corpak VTR				
	300	12.00	1.20	MS3030-P
Cosmed				
	Pony	4.80	1.30	MSL-3624
Covidien				
	Glide Scope	0.00	0.00	OEM ONLY
Criticare				
	100	12.00	2.60	L12V2.6F1-02
	501	4.80	0.60	MSL-5250
	503	6.00	1.30	L06V1.3F1-09
	503 Oximeter	6.00	1.20	MSL-5746
	503S Oximeter	6.00	1.20	MSL-5746
	504	6.00	5.00	L06V4.5F1-01
	504 DXP	6.00	3.20	L06V3.2F1-06
	5040 Pony	6.00	1.30	L06V1.3F1-09
	506	12.00	3.40	L12V3.4F1-04
	506 DXNP	6.00	5.00	L06V4.5F1-01
	506 NT3	6.00	7.20	L06V7.2F1-02
	507	12.00	3.40	L12V3.4F1-04
	507S, 508	12.00	2.80	L12V2.4F1-06
	602-14	7.20	2.70	MSL-5957
	8100 (newer)	12.00	2.10	L12V2.3C1-17
	8100 (older)	12.00	2.10	L12V2.3C1-17
	SPOT	6.00	1.30	L06V1.3F1-09
Critikon				
	1846	12.00	2.80	L12V2.4F1-06

COMPANY	MODEL	VOLTS	AH	BATTERY
Critikon				
	1846 Clock	3.60	0.90	MSL-5466
	8100	12.00	2.80	L12V2.4F1-06
	8700	12.00	2.10	L12V2.3C1-17
	9300	12.00	2.80	L12V2.4F1-06
	Compact T	12.00	2.10	L12V2.3C1-17
	Dinamap 8100	12.00	2.80	L12V2.4F1-06
	Dynamap mem	3.00	1.20	CR1/2AA-AX
	Dynamap Plus *	12.00	2.10	L12V2.3C1-17
	Dynamap Pro1000	0.00	0.00	OEM ONLY
	MPS	12.00	4.00	MSL-633153
	PRO 100	12.00	2.10	L12V2.3C1-17
	PRO300	12.00	2.10	L12V2.3C1-17
	PRO400	12.00	2.10	L12V2.3C1-17
	ProCare 300	6.00	3.00	L06V3.0F1-02
	ProCare 400	6.00	3.00	L06V3.0F1-02
Cross Vent				
	3	0.00	0.00	MSL-CV3
Cyber Power				
	685AVR	12.00	9.00	L12V9.0F2-04
Cyprus Lifts				
		12.00	5.00	L12V5.0F2-03
	1500	12.00	5.00	L12V5.0F2-03
Datascope				
	800 Monitors	12.00	2.10	L12V2.3C1-17
	Accutor	3.60	0.60	MSL-5488/T
	Accutor 4	12.00	2.20	MSL-5298
	Accutor III	12.00	2.60	L12V2.6F1-02
	Accutor IV	12.00	2.60	L12V2.6F1-02
	Accutor Plus	12.00	2.10	L12V2.3C1-17
	Accutor Plus	11.10	4.80	MSL-PM8000E
	Def DPD MRL	12.00	1.30	MSL-5715
	MD3	24.00	1.30	M300
	Passport *	12.00	2.10	L12V2.3C1-17
	System 97	12.00	18.00	L12V18B1-14
Datex-Ohmeda				
	017006	6.00	1.20	L06V1.2F1-09
	0690-1000-310	7.20	0.15	M9
	197230	2.80	0.12	MSL-6315
	3700	8.00	2.50	GL825
	3740	12.00	2.60	L12V2.6F1-02
	3770	9.60	1.40	MSL-3775-I
	3775	9.60	1.40	MSL-3775-I
	3800	8.00	2.70	L08V2.5F1-04
	3800/P	8.00	3.20	MSL-5921
	7100	6.00	5.00	L06V4.5F1-01
	7800 Vent	0.00	0.00	OEM ONLY
	7900	12.00	2.20	MSL-5899
	ADU *	12.00	2.60	L12V2.6F1-02

* May require multiple batteries

COMPANY	MODEL	VOLTS	AH	BATTERY
Datex-Ohmeda				
	Aestiva	12.00	2.20	MSL-5899
	Aestiva	12.00	2.60	L12V2.6F1-02
	AS3 Compact	24.00	1.70	MSL-5833NC
	AS3 PCB battery	0.00	0.00	OEM ONLY
	AS3-LITE	0.00	0.00	OEM ONLY
	AS5 B/U	0.00	0.00	OEM ONLY
	Cardio-Cap	3.00	1.20	CR1/2AA-AX
	CardioCap 5	12.00	2.60	L12V2.6F1-04
	Cardiocap5 B/U	0.00	0.00	OEM ONLY
	EE090178	24.00	4.00	MSL-5010
	FCM Rec100 Monitor	24.00	4.00	MSL-5010
	Giraffe	8.40	0.10	9V-NIMH
	Incubator	12.00	26.00	L12V26B1-16
	M4T28-BR12SHI	2.80	0.12	MSL-6315
	Modulus II	4.80	1.30	MSL-5100
	Modulus II+	6.00	12.00	L06V12F1-02
	O2 MRI	0.00	0.00	OEM ONLY
	Panda Baby Warmer	8.40	0.10	9V-NIMH
	RGM5250	3.00	0.00	CR2450
	S3 (Set)	24.00	4.00	MSL-5010
	S5 (Set)	24.00	4.00	MSL-5010
Delfi (Delphi)				
	PTS2	0.00	0.00	OEM ONLY
Delta Intellisyste				
	Alarm Panel	3.60	0.00	MSL-2216
Devilbiss				
	7304A, D	12.00	1.30	L12V1.3F1-04
	7304C	12.00	1.30	L12V1.3F1-04
	7305P	12.00	2.50	MSL-5975
	7310PR-D	0.00	0.00	OEM ONLY
Diagnostic Ultra				
	BVI2500	6.00	5.00	0809-0012
	BVI5000	6.00	5.00	0809-0012
Dolphin				
	Lift	12.00	7.20	L12V7.2F2-04
Drager				
	2HR-4UC-MC-2	2.40	0.70	MSL-JM103
	5732354	0.00	0.00	OEM ONLY
	Baby Log VN500	24.00	2.40	MSL-8415290/R
	Caleo Incubator	3.00	1.50	CR12600SE
	Delta	14.80	4.50	MSL-SC6802
	Eveta 2, 4, XL	3.00	1.50	CR12600SE
	Evita 4 (Internal) *	12.00	3.40	L12V3.4F1-04
	Evita XI (External)	12.00	17.00	LC-RD1217P
	Fabius GS	12.00	3.40	L12V3.4F1-04
	Gamma X	0.00	0.00	OEM ONLY
	Infinity Delta XL	0.00	0.00	OEM ONLY
	Jaundice Meter JM103	2.40	0.70	MSL-JM103

COMPANY	MODEL	VOLTS	AH	BATTERY
Drager				
	JM105	0.00	0.00	OEM ONLY
	Julian *	12.00	7.20	L12V7.2F2-04
	M540 (MS25794)	0.00	0.00	OEM ONLY
	M540 Infinity	0.00	0.00	OEM ONLY
	MS14490	0.00	0.00	OEM ONLY
	MS18340	0.00	0.00	OEM ONLY
	MS29574	0.00	0.00	OEM ONLY
	MU24842	0.00	0.00	OEM ONLY
	NarcomedMRJ	3.00		CR2335/P
	TI-500	12.00	26.00	L12V26B1-16
	Vapor 2000	0.00	0.00	OEM ONLY
Duracell				
		6.00	0.00	7K67
Dynasty				
	JC1222	12.00	2.60	L12V2.6F1-02
Enraf-Nonius				
	463, 491 SonoPuls490	12.00	2.60	L12V2.6F1-02
Entech				
	Enteral Pump	4.00	3.00	L04V3.0F1-16
Ergolet				
	600	12.00	7.20	L12V7.2F1-04
Filac				
	1010, 528	8.40	0.20	MSL-5595
First Temp				
		2.40	2.20	MS222/T
	2000A	2.40	4.00	MS240/T
Fluke				
	BP7217	1.20	1.60	MSL-P160AS
Foxborough				
	Miran 2	7.20	0.00	MSL-FOX7.2
Fresenius				
	5008 *	12.00	7.20	L12V7.2F2-04
	Agilia MC	0.00	0.00	OEM ONLY
Gambro				
	Artis	12.00	7.20	L12V7.2F2-04
	Phoenix (6973440)	12.00	7.20	L12V7.2F2-04
GE/Marquette				
	2016989-002	0.00	0.00	OEM ONLY
	2017857-002	11.10	4.20	MSL-5071
	2023227-001	0.00	0.00	OEM ONLY
	2023582-029	0.00	0.00	OEM ONLY
	2031069-003	0.00	0.00	OEM ONLY
	2035701-001	0.00	0.00	OEM ONLY
	2036984-001	0.00	0.00	OEM ONLY
	2036984-001	0.00	0.00	OEM ONLY
	2062895-001	0.00	0.00	OEM ONLY
	2069650-001	3.00	0.10	CR2032-PA
	900 Defib Recorder	24.00	1.20	MSL-5278

* May require multiple batteries

COMPANY	MODEL	VOLTS	AH	BATTERY
GE/Marquette	Aisy Anesthesia	12.00	12.00	L12V12F2-04
	AMX-IV G26EP (Req 9)	12.00	26.00	MSL-5970
	AMX-IV Pana. (Req 9)	12.00	28.00	L12V28B2-15
	AMX-IV Peak (Req 9)	12.00	26.00	L12V26B1-16
	Carescape B450	0.00	0.00	OEM ONLY
	Carescape B650	0.00	0.00	OEM ONLY
	Carescape V100	6.00	3.00	L06V3.0F1-02
	Carescapre PDM	0.00	0.00	OEM ONLY
	C-ARM OC9900	0.00	0.00	OEM ONLY
	Dash 2500	0.00	0.00	OEM ONLY
	DASH 3000 *	11.10	4.20	MSL-5071
	Dash 4000	11.10	4.20	MSL-5071
	Dash 5000	11.10	4.20	MSL-5071
	Defib 46	24.00	1.30	MSL-5014
	Eagle	0.00	0.00	OEM ONLY
	FLEX-3S2P	0.00	0.00	OEM ONLY
	FLEX-3S3P	0.00	0.00	OEM ONLY
	Logiq E US Notebook	0.00	0.00	OEM ONLY
	M1168350	3.00		CR1632
	MAC 12	3.60	0.60	MSL-5653
	MAC 1600	0.00	0.00	OEM ONLY
	MAC 5000	0.00	0.00	OEM ONLY
	MAC 5500	0.00	0.00	OEM ONLY
	MAC 6	0.00	0.00	OEM ONLY
	MAC PC	0.00	0.00	OEM ONLY
	Mac Vu	0.00	0.00	OEM ONLY
	MAC1200	18.00	1.30	MSL-5877
	Mac2 7100 *	3.60	0.15	MEMPAK3.6V
	Mac-8	0.00	0.00	OEM ONLY
	PDM Mini Dash	0.00	0.00	OEM ONLY
	ProCare 400	6.00	3.40	L06V3.4F1-01
	Tram	0.00	2.80	403258-001
	TRAM PACII	12.00	2.40	MSL-5494
	TRAM7200	12.00	2.00	L12V2.0C1-17
Gomco	178 Aspirat	12.00	4.50	L12V5.0F1-03
	G180 Suction	12.00	4.50	L12V5.0F1-03
Gould	SP2204 Recorder	12.00	5.00	GL0800-0109
Graseby	3000	12.00	1.30	L12V1.3F1-04
	3100	2.00	2.50	3XGL0810-0004
	3300	2.00	2.50	3XGL0810-0004
	3400	2.00	2.50	3XGL0810-0004
Grieshaber	Cutter	7.20	0.60	MSL-5760
Guldmann	Lift	12.00	2.70	L12V2.7F1-03

* May require multiple batteries

COMPANY	MODEL	VOLTS	AH	BATTERY
Hamilton (vents)				
	Galileo	12.00	4.50	L12V5.0F1-03
Handicare				
	C1000 *	12.00	4.50	L12V5.0F1-03
	C-450	12.00	4.50	L12V5.0F1-03
	C-625	12.00	4.50	L12V5.0F1-03
Hawker/Gates				
		8.00	5.00	0800-0027
	8416	6.00	8.00	MSL-NPB4000
	D Size	2.00	2.50	GL0810-0004
	X Size	2.00	5.00	GL0800-0004
Healthdyne				
	Smart Monitor 970S	6.00	2.30	MSL-5492
HealthOmeter				
	2200	0.00	0.00	OEM ONLY
	2600KLS	0.00	0.00	OEM ONLY
	451KL	6.00	3.00	L06V3.0F1-02
	551KL	6.00	2.80	L06V2.4F1-06
	552KL	6.00	3.00	L06V3.0F1-02
	592KL	6.00	3.20	L06V3.2F1-06
	595KL	6.00	3.20	L06V3.2F1-06
	63719	6.00	3.00	L06V3.0F1-02
	970S	6.00	2.30	MSL-5492
Heine				
	Accubox II	0.00	0.00	OEM ONLY
	Beta-200	3.50	0.75	MSL-BETA-R
	F22-804	2.40	0.70	MSL-F22-804
	X-02.99.314	0.00	0.00	OEM ONLY
	X-02.99.830	3.50	0.75	X-02.99.382
	X-02.99.832	3.50	0.75	X-02.99.382
	X-04.99.623	0.00	0.00	OEM ONLY
	X-04.99.624	0.00	0.00	OEM ONLY
Heska				
	4400	6.00	0.80	MSL-4400
	4402L	0.00	0.00	MSL-4402L
Hilger				
	2R/H2	8.40	2.30	MSL-2R
	Nerve St H2	15.00	0.00	MSL-3159
Hillrom				
	Progressa Bed*	12.00	5.00	L12V5.0F2-03
	VersaCare	12.00	7.20	L12V7.2F2-04
Hokanson				
	MD6VR	14.40	0.70	MSL-MD6R
Hospira/Abbott				
	0100-0304	4.80	2.10	MSL-0100-0304
	Acclaim	8.00	2.70	L08V2.5F1-04
	Acclaim	8.00	2.70	L08V2.5F1-04
	Flexi Companion	4.00	1.00	L04V1.0C1-17
	Flexi Flo	6.00	2.50	MSL-5223

* May require multiple batteries

COMPANY	MODEL	VOLTS	AH	BATTERY
Hospira/Abbott				
	Flexi Patrol	6.00	5.00	L06V4.5F1-01
	Flexiflo III	6.00	2.80	L06V2.4F1-06
	Istat Printer	4.80	1.60	MSL-6192
	Lifecare 5000	8.00	2.70	L08V2.5F1-04
	Lifecare III	6.00	2.50	GL4018
	Lifecare IV	8.00	2.70	L08V2.5F1-04
	Lifecare PCA	8.00	2.70	L08V2.5F1-04
	Lifecare PCA	8.00	2.70	L08V2.5F1-04
	LIN337-001	0.00	0.00	OEM ONLY
	MCP9819-065	4.80	1.60	MSL-6192
	MJ09.01	4.80	1.60	MSL-6192
	Patrol EXT Pump	0.00	0.00	OEM ONLY
	PCA Plus4100	8.00	2.70	L08V2.5F1-04
	PCA Printer	4.80	0.60	MSL-AA4.8
	Plum	8.00	2.70	L08V2.5F1-04
	Plum (Plug)	8.00	2.70	L08V2.5F1-04
	Plum 360	6.00	5.00	L06V4.5F1-01
	Plum A Plus	6.00	5.00	L06V4.5F1-01
	Plum XL	8.00	2.50	MSL-5875
	Sapphire (16354-01)	0.00	0.00	OEM ONLY
	Simbiq	0.00	0.00	OEM ONLY
Housted				
	Elect Bed	12.00	1.20	L12V1.2F1-04
Hoyer				
	HPL-P	12.00	2.70	L12V2.7F1-03
Human Scale				
	T5	12.00	12.00	L12V12F2-04
IMED				
	900 Series	6.00	12.00	L06V12F1-02
	Gemini PC1	6.00	7.20	L06V7.2F1-02
	Gemini PC2	6.00	12.00	L06V12F1-02
	Gemini PC4	12.00	3.40	L12V3.4F1-04
	Gemini TX *	6.00	3.40	L06V3.4F1-01
	PC Series mem b/u	4.80	0.70	MSL-I4.8
Imex				
	CT+	12.00	0.25	MSL-5674NMH
	Doppler CtPI	0.00	0.00	OEM ONLY
Impact				
	302 Suction	6.00	3.00	L06V3.0F1-02
	305 Suction *	6.00	5.00	L06V4.5F1-01
	320 Resp.	6.00	1.30	L06V1.3F1-09
	321NOH	12.00	2.00	L12V2.0F1-08
	750 Vent	12.00	3.20	MSL-5642
	VAC-PAC	6.00	3.20	L06V3.2F1-06
Infasonics				
	Adult Star	12.00	7.20	L12V7.2F1-04
	Infant Star	10.00	2.50	MSL-5342
Innova				

COMPANY	MODEL	VOLTS	AH	BATTERY
Innova				
	OCU-CEL XL	3.00	0.00	DL1/3N
	Tonopen XL	3.00	0.00	DL1/3N
Instromedix				
	Express 3	3.00	0.35	LTC-3PN
Integra / Gambro				
	Dialysis	0.00	0.00	OEM ONLY
	Phoenix	0.00	0.00	OEM ONLY
Intermec				
		1.00	1.00	N1000AAE
International Bio				
	Voyager 626	12.00	26.00	MSL-5970
Intl Technodyne				
	801	7.20	2.30	MSL-5470
	Hemochron 8430	7.20	5.00	MSL-5083
	Hemochron Signature Jr	8.40	1.65	MSL-5974NIMH
Invacare				
	Reliant 450	12.00	2.70	L12V2.7F1-03
	RPL4501	12.00	2.70	L12V2.7F1-03
	RSA660	12.00	2.70	L12V2.7F1-03
Invensys				
	DMS-BATT-1	3.60	0.10	MS72500
Invivo				
	4500	8.00	2.50	90-M8V4500
	Omega	12.00	2.60	L12V2.6F1-02
	Scout	4.80	1.50	MSL-5711
Jaeger				
	Spiroprot	0.00	0.00	OEM ONLY
Jones				
	Spirometer	7.20	0.60	MSL-5332
Jostra				
	HL-20	12.00	30.00	A512/30G6
	HL-30	12.00	28.00	LC-X1228AP
	RFC20-970	0.00	0.00	OEM ONLY
KCI				
	414003	0.00	0.00	OEM ONLY
	Activac	0.00	0.00	OEM ONLY
	Infovac	0.00	0.00	OEM ONLY
Krups				
	Jetpack	2.40	0.60	MSL-2/3Ax2
Laerdal				
	3000	12.00	2.00	L12V2.0C1-17
	7050 Suction	12.00	1.20	STR5
	7800	12.00	2.10	L12V2.3C1-17
	790010	12.00	1.20	STR5
	884301	12.00	2.50	MSL-5975
	886113	0.00	0.00	OEM ONLY
	Compac Suctn	12.00	1.30	L12V1.3F1-04
	HS 3000	12.00	2.00	L12V2.0C1-17

COMPANY	MODEL	VOLTS	AH	BATTERY
Laerdal				
	Impact	12.00	5.00	MSL-5269
	Impact 305 *	6.00	5.00	L06V4.5F1-01
	LCSU (880020)	12.00	2.50	MSL-5975
	LCSU 2	12.00	2.50	MSL-5975
	LCSU3	0.00	0.00	OEM ONLY
	LCSU4	0.00	0.00	OEM ONLY
	LSU Suction	12.00	2.10	L12V2.3C1-17
	Memory BU	3.00	0.60	MSL-5851
	Suction	12.00	1.20	L12V1.2F1-04
Leica				
	Footswitch (ER34615)	3.60		LS33600
Lifecare				
	PLV100	12.00	2.50	MSL-35-277
Lifeline				
	H101	6.00	1.00	L06V1.0F1-01
Liko (Patient Lift)				
	242	12.00	2.60	L12V2.6F1-04
	Likorall R2R *	12.00	2.60	L12V2.6F1-04
	Sabina II*	12.00	2.70	L12V2.7F1-03
Lionville				
	Medcart 34015-2	6.00	7.20	L06V7.2F1-02
	Medcart 400	6.00	7.20	L06V7.2F1-02
	Medcart 800	6.00	20.00	MSL-6047
	Medcart i.Point3	12.00	18.00	L12V18B1-14
Lithium				
		3.00	0.00	CR2450
		3.00	1.20	CR1/2AA
		3.60		LS14500BA
		9.00	0.00	CR9V
		3.00	1.80	CR17335 SE/T2
		3.00	0.10	CR2325
		3.00	1.50	CR12600SE
		3.00	0.15	CR2025
		3.00	1.40	CR2/3A
		3.00	0.10	CR2032-PA
		3.00	1.20	CR1/2AA-AX
	1/2 AA with Tabs	3.60	1.20	LS14250/T
	AA	1.50	2.70	CRAA-1.5
Maquet				
	OR Table *	12.00	12.00	L12V12F2-04
	Servo I (6487180)	0.00	0.00	OEM ONLY
Masimo				
	Rad 87	6.00	1.30	L06V1.3F1-09
	Rad8	6.00	1.30	L06V1.3F1-09
	Radical 7 Handheld	0.00	0.00	OEM ONLY
Matrix				
	Life Def+	12.00	2.30	MSL-5701 (REF)
MDE				

* May require multiple batteries

COMPANY	MODEL	VOLTS	AH	BATTERY
MDE				
	100, 200	12.00	2.00	L12V2.0C1-17
	100T	12.00	3.20	MSL-5441
	300A	12.00	2.10	L12V2.3C1-17
	300A 300E	12.00	2.10	L12V2.3C1-17
	B/U	0.00	0.00	OEM ONLY
	Escort II	12.00	2.10	L12V2.3C1-17
	Prism	12.00	7.00	MSL-5045
Medasonics				
	D8, D9, D10	0.00	0.00	OEM ONLY
Medela				
	600.0806	12.00	2.10	MSL-6384
	919.0012	0.00	0.00	OEM ONLY
	Vario 18	0.00	0.00	OEM ONLY
	Vario 8	0.00	0.00	OEM ONLY
Medfusion Syste				
	2001	7.20	1.50	MSL-5410
	2010 Syringe	7.20	1.50	MSL-5410
	22-4050-51	0.00	0.00	OEM ONLY
	3010	7.20	2.10	MSL-5181
	3500 (w/PCB)	7.20	2.10	MSL-5181
	4000	0.00	0.00	OEM ONLY
	7000	7.20	4.00	MSL-7000
	III	3.60	1.80	MSL-5624-I
	Protege	7.20	2.10	MSL-5181
	Protégé 3010	7.20	2.10	MSL-5181
Medi Lifter				
	3	12.00	5.00	L12V5.0F2-03
Medical Positioni				
	Echo Bed	6.00	1.20	L06V1.2F1-09
Medi-Man				
	300*	12.00	4.50	L12V5.0F1-03
	4000 *	12.00	5.00	L12V5.0F2-03
	4050	12.00	5.00	L12V5.0F2-03
	5000 *	12.00	7.20	L12V7.2F1-04
	7000 *	12.00	5.00	L12V5.0F2-03
	RDP-3*	12.00	4.50	L12V5.0F1-03
Micro Medical				
	A200 Spiro	7.20	0.60	MSL-5085
	MLA5009	7.20	0.60	MSL-6144
Midmark				
	7100 *	12.00	12.00	L12V12F1-04
Mindray				
	PM-8000	0.00	0.00	OEM ONLY
	PM-9000	0.00	0.00	OEM ONLY
	VS800	0.00	0.00	OEM ONLY
	Zoneare Pro	0.00	0.00	OEM ONLY
Minogue Medica				
	E09711	12.00	2.00	L12V2.0C2-17

* May require multiple batteries

COMPANY	MODEL	VOLTS	AH	BATTERY
Minogue Medica	GPS Navigator	12.00	2.00	L12V2.0C2-17
Minolta	JM101	4.80	0.20	MS101
	JM102	0.00	0.00	OEM ONLY
	JM103	2.40	0.70	MSL-JM103
	SAT	0.00	0.00	MS08
	Sat	7.20	1.00	90-M10265
Mira	Mira TR3000	24.00	1.20	MB5615
	TR4000, 3000	0.00	0.00	OEM ONLY
Monitral	Dialysis	12.00	0.70	MSL-9002
Mortara Instrum	4800-011	0.00	0.00	OEM ONLY
	4800-017	0.00	0.00	OEM ONLY
	ELI 10	0.00	0.00	OEM ONLY
	ELI 100	6.00	1.00	L06V1.0F1-01
	ELI 280, 350	12.00	9.00	L12V9.0F1-04
	ELI 350 (4800-010)	12.00	9.00	L12V9.0F2-04
	ELI 380	0.00	0.00	OEM ONLY
MRL	90 PortapakD	12.00	1.50	MSL-MRL-90
	9008 Def	12.00	1.50	MSL-5528
	Porta Pak 80	12.00	2.50	MSL-MRL-80
	Porta Pak 90	12.00	2.20	MSL-5229
	SmartPac 360	12.00	1.30	MSL-5554
NA Draeger	Oxilog 2000	7.20	1.20	MSL-5938
Narco	Scale	12.00	2.60	L12V2.6F1-02
	Scale N10	12.00	2.60	L12V2.6F1-02
Nashua	Med Cart	12.00	4.50	L12V5.0F1-03
Nellcor	N10	6.00	0.50	MS71764
	N100	12.00	2.50	GL0810-0115
	N1000	6.00	5.00	L06V4.5F1-01
	N180	12.00	2.10	L12V2.3C1-17
	N200	6.00	2.80	L06V2.4F1-06
	N290	6.00	5.00	L06V4.5F1-01
	N3000	12.00	2.00	L12V2.0F1-08
	N-3000 All	12.00	2.00	L12V2.0F1-08
	N3100	3.00	0.00	CR1225
	N390	6.00	5.00	L06V4.5F1-01
	N550	12.00	1.20	L12V1.2F1-04
	N-560	9.60	3.80	MSL-6093
	N595	6.00	5.00	L06V4.5F1-01
	N-600	6.00	5.00	L06V4.5F1-01

COMPANY	MODEL	VOLTS	AH	BATTERY
Nellcor	N-85	7.20	2.70	MSL-5914-I
	NPB	6.00	5.00	L06V4.5F1-01
	NPB190	6.00	5.00	L06V4.5F1-01
	NPB295/395/39xx	6.00	5.00	L06V4.5F1-01
	NPB4000	6.00	8.00	MSL-NPB4000
Nicad		1.20	0.00	MSL-4/5A/T
		2.40	0.00	MSL-2/3A/TX2
		3.60	0.00	MSL-2/3AA/TX3
	1/2 AA	1.20	0.20	N270AA
	1/2 C	1.20	0.00	MSL-1/2 C/T
	1/2 Sub C	1.20	0.00	MSL-1/2SC/T
	3 AA Tabs	0.00	0.00	MSL-AA/TX3
	4 C w Tabs	0.00	0.00	MSL-CX4/T
	5 AA Tabs	0.00	0.00	MSL-AA/TX5
	7 C w Tabs	0.00	0.00	MSL-CX7/T
	8 AA Tabs	0.00	0.00	MSL-AA/TX8
	AA	1.20	0.00	MSL-AA
	AA Tabs	1.20	0.00	MSL-AA/T
	AAA	1.20	0.00	N200AAA
	AAA Tabs	1.20	0.00	MSL-AAA/T
	C	1.20	1.70	MSL-C1700
	D	1.20	4.00	MSL-D
	F	0.00	0.00	MSL-F
	Full AA	1.20	0.00	MSL-AA
	Full C	1.20	2.80	MSL-C2800
	Full D	1.20	4.30	MSL-D4300
	Sub C	1.20	1.80	MSL-SC
	Sub C	0.00	0.00	MSL-SC4/5/T
	Sub C	1.20	0.00	MSL-SC1400
	Sub C	1.20	0.00	MSL-SC1800
	Sub C w Tabs	1.20	1.80	MSL-SC/T
	Sub C w Tabs	0.00	0.00	MSL-SC1700T
Nihon Kohden	2101	6.00	3.40	L06V3.4F1-01
	6201A	12.00	2.00	L12V2.0C1-17
	Cardio Fax	3.60		LS14250BA
	Defib 101	12.00	1.70	NK101-DEF
	Lifescop N	9.60	2.10	MSL-6042
	TEC 7200	12.00	2.10	L12V2.3C1-17
NiMH		3.60	0.60	3/V70H-PC
	7.2 Volt	7.20	0.15	M9
	9 Volt	8.40	0.10	9V-NIMH
	AA	1.20	0.00	MSL-AANIMH
	AA Tabs	1.20	0.00	MSL-AANIMH/T
	AA Tabs	1.20	0.00	MSL-AA1500/T
	AA/Tabs	1.20	0.00	MSL-AA1000/T

COMPANY	MODEL	VOLTS	AH	BATTERY
Nonin	2120	7.20	3.80	MSL-2120
	7500	0.00	0.00	OEM ONLY
	8000/8800	6.00	1.50	MSL-5353/T
	8600	6.00	1.50	MSL-5353/T
	9600	7.20	3.80	MSL-2120
Novamatrix	1265/515	12.00	2.60	L12V2.6F1-02
	500	6.00	3.20	L06V3.2F1-03
	520A	12.00	2.80	L12V2.4F1-06
	8200	6.00	3.20	L06V3.2F1-06
	Tidalwave	7.20	3.80	DR30AA
Ohio MED6008	AirVac	12.00	26.00	L12V26B1-16
Olympic	56320	7.20	1.80	MSL-5830
	57193	6.00	1.50	MSL-5296-I
	Smart Scl	7.20	1.80	MSL-5830
Olympus	CLV-U10	3.00	1.80	MSL-5779
	CLV-U20	3.60	2.10	MSL-5084B
	CV-180	0.00	0.00	OEM ONLY
	M4T28-BR12SH	0.00	0.00	OEM ONLY
OmniCell	ER34615	3.60		LS33600
Omron Healthcar	48H907N	4.80	1.60	MSL-6209
	HEM-907XL	4.80	1.60	MSL-6209
Palco	Sat 340	0.00	0.00	OEM ONLY
Parks	811 Doppler	12.00	0.20	M811
	811B	12.00	1.20	L12V1.2F1-04
	915L	12.00	1.20	L12V1.2F1-04
Patient Tech	Survalent	6.00	0.20	MSL-5412
Penlon	AM-1000	12.00	2.60	L12V2.6F1-02
	AV 500	3.60	0.10	MSL-5552
	MRI Laryngoscope	0.00	0.00	OEM ONLY
	MRI Laryngoscope	0.00	0.00	OEM ONLY
Philips/HP	1700 ECG	18.00	7.20	MSL-1721A
	200i EKG	0.00	0.00	OEM ONLY
	40488A	12.00	2.10	L12V2.3C1-17
	43110	12.00	2.50	GL43100
	453561731923	0.00	0.00	OEM ONLY
	4745A	3.60	0.60	MSL-5180
	78353	24.00	2.20	MS78353

COMPANY	MODEL	VOLTS	AH	BATTERY
Philips/HP				
	78554	7.20	0.60	MS78554
	78560	12.00	2.20	MS78560
	78619 78620	12.00	4.50	MS78619
	78620	12.00	4.00	MS78620
	78650	3.60	0.60	MSL-5531I
	78660	12.00	2.20	MS78660
	78670	0.00	0.00	OEM ONLY
	8026A	6.00	0.20	MS5R225
	8040	0.00	0.00	OEM ONLY
	989803129131	0.00	0.00	OEM ONLY
	989803130151	0.00	0.00	OEM ONLY
	989803135861	10.80	6.00	MSL-6139
	98980315161	0.00	0.00	OEM ONLY
	989803166291	0.00	0.00	OEM ONLY
	989803170371	0.00	0.00	OEM ONLY
	989803174131	0.00	0.00	OEM ONLY
	989803194541	11.10	7.80	MSL-6181
	A3	6.00	8.00	MSL-NPB4000
	C1	6.00	5.00	L06V4.5F1-01
	C3	12.00	4.50	L12V5.0F1-03
	CdmstrM1722A	12.00	4.00	MSL-1758A-INT
	Code Master	6.00	5.00	L06V4.5F1-01
	Codemstr XL	6.00	5.00	L06V4.5F1-01
	CX50	0.00	0.00	OEM ONLY
	Data Mgmt	7.20	0.70	MSL-78554-I
	FM20, FM30	0.00	0.00	OEM ONLY
	FR3	0.00	0.00	OEM ONLY
	Heartstart M3516A	0.00	0.00	OEM ONLY
	M1125A	0.00	0.00	OEM ONLY
	M1205A	12.00	2.10	L12V2.3C1-17
	M1275A *	12.00	2.10	L12V2.3C1-17
	M1276A	12.00	2.10	L12V2.3C1-17
	M13560, clock *	3.00	0.70	MSL-5709
	M1770	0.00	0.00	OEM ONLY
	M2460A	0.00	0.00	OEM ONLY
	M2476B	0.00	0.00	OEM ONLY
	M2477B	12.00	4.00	MSL-M2477B
	M2636 Telemon	0.00	0.00	OEM ONLY
	M2636B	0.00	0.00	OEM ONLY
	M2720-64101	0.00	0.00	OEM ONLY
	M3	12.00	3.50	MSL-M3046A
	M3002A	0.00	0.00	OEM ONLY
	M3046A	12.00	3.50	MSL-M3046A
	M3516A	0.00	0.00	OEM ONLY
	M3535/M3536	0.00	0.00	OEM ONLY
	M3538A	0.00	0.00	OEM ONLY
	M4605A	10.80	6.00	MSL-6139
	M4607A	0.00	0.00	OEM ONLY

* May require multiple batteries

COMPANY	MODEL	VOLTS	AH	BATTERY	
Philips/HP	M4735A	0.00	0.00	OEM ONLY	
	M4790A	10.80	3.80	MSL-5981	
	Merlin	12.00	2.10	L12V2.3C1-17	
	MP2	10.80	1.00	MSL-4607A/R	
	MP20	10.80	6.00	MSL-6139	
	MP30	10.80	6.00	MSL-6139	
	MRX	0.00	0.00	OEM ONLY	
	MX40	0.00	0.00	OEM ONLY	
	N1030HP	0.00	0.00	OEM ONLY	
	Omni-Care Tr	12.00	2.10	L12V2.3C1-17	
	Page Writer	6.00	7.20	L06V7.2F1-02	
	Pagewriter Touch	0.00	0.00	OEM ONLY	
	Pgwtr 1700	18.00	7.20	MSL-1721A	
	Sparq	0.00	0.00	OEM ONLY	
	SureSigns VS2+	0.00	0.00	OEM ONLY	
	TC70	0.00	0.00	OEM ONLY	
	Telemon M2636A/B/C	10.80	3.80	MSL-5981	
	TRIM3	0.00	0.00	OEM ONLY	
	V26	12.00	3.50	MSL-M3046A	
	V60	0.00	0.00	OEM ONLY	
	VM3, VM4, VM6	11.10	7.80	MSL-6181	
	VS1	6.00	5.00	MSL-6014	
	VS2	11.10	7.80	MSL-6181	
	VS2+	0.00	0.00	OEM ONLY	
	VS3	11.10	7.80	MSL-6181	
	VS4	11.10	7.80	MSL-6181	
	VS6	11.10	7.80	MSL-6181	
	XLI	18.00	7.20	MSL-1721A	
	Physio Control	11141-000149	0.00	0.00	OEM ONLY
		300	12.00	3.40	LP300
		3200497-000	12.00	3.00	MSL-5042
		First Medic	12.00	2.00	L12V2.0C2-17
LifeStat 100		8.00	2.50	GL825-100-2P	
LifeStat 1600		8.00	2.50	GL825-200	
LP-10		12.00	1.00	LP5-FP	
LP1000		0.00	0.00	OEM ONLY	
LP12		0.00	0.00	OEM ONLY	
LP12-SLA		0.00	0.00	OEM ONLY	
LP15		0.00	0.00	OEM ONLY	
LP20		0.00	0.00	OEM ONLY	
LP20 Memory		3.00	0.10	CR2032-PA	
LP200		16.00	2.50	MSL-5301	
LP20E		0.00	0.00	OEM ONLY	
LP4		12.00	4.30	LP4	
LP5		12.00	1.00	LP5	
LP500		0.00	0.00	OEM ONLY	
LP6 Defib		12.00	1.00	LP6/200	

COMPANY	MODEL	VOLTS	AH	BATTERY
Physio Control				
	LP6 Monitor	14.00	1.50	LP6/100
	LP7 Defib	12.00	1.00	LP6/200
	LP7 Monitor	14.00	1.50	LP6/100
	LP8 Defib	0.00	0.00	OEM ONLY
	LP8 Mon	0.00	0.00	OEM ONLY
	LP9	16.00	2.80	LP9-S
	VSM2 Monitor	14.40	2.30	MSL-5029
	VSM5 Clock	3.60	0.10	MSL-5333
Pipett				
	Pipett	6.00	0.15	5V150H
Plantronics				
	65358-01	3.70	0.23	MSL-BATT-CS50
Power Battery C				
	PRC-1230X	12.00	30.00	L12V35B1-15
Precision Medica				
	EasyGo	12.00	2.70	L12V2.7F1-03
	PM65	12.00	2.70	L12V2.7F1-03
Prism Medical				
	F-550	12.00	2.70	L12V2.7F1-03
	Lift	12.00	2.60	L12V2.6F1-02
	Lift	12.00	2.70	L12V2.7F1-03
	P440*	12.00	2.10	L12V2.3C1-17
	Transactive Plus	12.00	7.20	L12V7.2F1-04
Protocol				
		8.00	2.70	L08V2.5F1-04
	ProPaq mem B/U	3.00	0.10	CR2032-PA
	ProPaq w Oximetry	8.00	5.40	MSL-PRO-S
Pulmonetics				
	LTV1000ext	12.00	12.00	L12V12F2-04
	LTV1000int	12.00	4.50	MSL-5117
	LTV1200 Internal	12.00	5.00	MSL-6212
Pulsair				
	OMS50	6.00	0.20	MSL-OMS50
Puritan Bennett				
	2801 Vent	12.00	4.00	MSL-5627
	7200 *	2.00	5.00	GL0800-0004
	840	24.00	7.00	MSL-5110
	LP10 Vent	12.00	2.80	L12V2.4F1-06
	PB100	3.60	0.60	MSL-5677
Radiometer				
	TCM3 Monitor	7.20	2.30	MSL-TCM-3
Rauland				
	BK2000	14.40	1.30	MSL-5705
	Nurse Call	0.00	0.00	OEM ONLY
Reichert				
	Avia Tono-Pen	0.00	0.00	OEM ONLY
	PN230577	0.00	0.00	OEM ONLY
	Tono-Pen XL *	3.00		CR1/3N

COMPANY	MODEL	VOLTS	AH	BATTERY
Respironics				
	1006005	2.80	0.12	MSL-6315
	1058036	0.00	0.00	OEM ONLY
	900-102	0.00	0.00	OEM ONLY
	Billichek	6.00	0.70	MSL-6090/R
	BiPap Vision	0.00	0.00	OEM ONLY
	EverGo	0.00	0.00	OEM ONLY
	ILS-100 / ILS-101	0.00	0.00	OEM ONLY
	Infinity	0.00	0.00	OEM ONLY
	Remstar Auto Bipap	3.00	0.10	CR2032-PA
	Smart Monitor	6.00	2.00	MSL-5996
Rice Lake				
	Scale 140-10-7	6.00	1.20	L06V1.2F1-09
RMD Instrument				
	Camera E097113	12.00	2.00	L12V2.0C2-17
	Camera E097113	12.00	2.00	L12V2.0C2-17
Roche				
	Inform	3.70	1.80	MSL-6298R
Rolm				
	NCS2000	12.00	1.70	MSL-NCS2000
Rubbermaid				
	9M09	0.00	0.00	OEM ONLY
	9M2109	0.00	0.00	OEM ONLY
Sanyo				
	5HR-AAC	6.00	1.60	MSL-5HR-AAC
	MDF-C8V	6.00	1.60	MSL-6445
	MDF-U32V	6.00	1.60	MSL-6445
	NIC0262	6.00	1.60	MSL-5HR-AAC
Scale-Tronix				
	20009	6.00	0.75	M2001
	6006	6.00	0.75	M2001
	Scale 4800	6.00	0.75	M2001
	Scale 5005	6.00	0.75	M2001
	Scale 6102	7.20	1.80	MSL-6006
Schiller				
	AT-10	12.00	0.80	L12V0.8W2-11
	AT-102	12.00	2.60	L12V2.6F1-02
	AT2, AT2 Plus	12.00	2.60	L12V2.6F1-02
	AT-5	9.60	1.60	MSL-5824
	ECG CS100	12.00	3.40	L12V3.4F1-04
	SP 10 Spiro *	12.00	0.80	L12V0.8W2-11
	SP200	12.00	2.60	L12V2.6F1-02
Scio-Teck				
	918850	12.00	3.40	L12V3.4F1-04
Seca				
	664	7.20	0.70	MSL-5784
	727	9.60	0.60	MSL-5750
	728	0.00	0.00	OEM ONLY
	944, 954	7.20	0.70	MSL-5784

COMPANY	MODEL	VOLTS	AH	BATTERY
Seca				
	CNH-130AAH-F2X3	0.00	0.00	OEM ONLY
Secrist				
	Millenium	12.00	2.10	MSL-6023
Sensormedics				
		9.60	2.20	MS765563
Shampaine				
	5100 *	12.00	26.00	L12V26B1-16
Sherwood				
	224 Kangaroo	7.20	1.50	MSL-5284
	324 Enhanced (010263	7.20	1.60	MSL-5997
	324 Kangaroo	7.20	1.50	MSL-5284
	520 / 2	8.40	0.20	MSL-5155
	Kangaroo 924	7.20	1.60	MSL-5070
	Kangaroo Control	7.20	1.60	MSL-5070
	Kangaroo E Pump	4.80	3.80	MSL-6094
	Pet Ambulat	3.60	1.80	MSL-5626
Shimatzu X-Ray				
	MC125L-50	12.00	30.00	L12V35B1-15
Siemens				
	300 Vent *	12.00	2.00	L12V2.0F1-09
	6002XL	12.00	2.00	L12V2.0C2-17
	64-87-180 (w/o Gauge	0.00	0.00	OEM ONLY
	6802 Excel	12.00	2.00	L12V2.0C2-17
	9000XL	0.00	0.00	OEM ONLY
	9000XL	12.00	2.10	L12V2.3C1-17
	BD30148	0.00	0.00	OEM ONLY
	Servo I	0.00	0.00	OEM ONLY
Sigma				
	6000	12.00	3.40	L12V3.4F1-04
	8000	12.00	2.50	MSL-5896
Site Rite				
	Lab 8001C	0.00	0.00	OEM ONLY
Skytron				
	1700 Essentia	12.00	18.00	L12V18B1-14
	3501-1	12.00	18.00	L12V18B1-14
	3501B	12.00	18.00	PS-12180F2
	6500 *	12.00	18.00	L12V18B1-14
	6700B Hercules	12.00	18.00	L12V18B1-14
	6701	12.00	18.00	L12V18B1-14
Smith Medical				
	Omnifuse	3.60	0.60	MSL-5488/T
Sonicaid				
	205	8.40	0.20	MSL-5233
	D205	9.60	0.50	MS8001
Sonosite				
	180 Plus	0.00	0.00	OEM ONLY
	MicroMax	0.00	0.00	OEM ONLY
	P07168-02	0.00	0.00	OEM ONLY

COMPANY	MODEL	VOLTS	AH	BATTERY
Sonosite				
	P0773	0.00	0.00	OEM ONLY
	S-ICU	0.00	0.00	OEM ONLY
Sony				
	Lithium	0.00	0.00	DR35S
Soren				
	S5	12.00	20.00	LC-X1220P
Spacelabs				
	1050 *	12.00	2.40	MSL-5389NIMH
	1050 Portable Monito	3.00	1.20	MSL-5062
	1050 Portable Monito	12.00	2.10	L12V2.3C1-17
	4045	24.00	1.30	MSL-5340
	600 Series	18.00	2.50	GL600-0016
	600A Series	18.00	3.20	L06V3.2F1-06X3
	90207	3.00	0.00	DL1/3N
	90305	8.00	3.50	L08V3.2F1-06
	90308 *	12.00	2.10	L12V2.3C1-17
	90309 *	12.00	2.10	L12V2.3C1-17
	90367 *	12.00	2.10	L12V2.3C1-17
	90479 Telemetry Rec.	12.00	3.20	MSL-90479
	90600 Series	18.00	3.20	L06V3.2F1-06X3
	90604A	6.00	3.20	L06V3.2F1-06
	90724-BLOCK	18.00	3.20	MSL-90724
	91387 Ultraview Mon.	3.00	1.20	MSL-5062
	91387 Ultraview Mon.	12.00	0.70	MSL-5943
	91390	0.00	0.00	OEM ONLY
	Medic 6			N/A
	MOM9400	12.00	2.10	L12V2.3C1-17
	PC Express	12.00	2.10	L12V2.3C1-17
	PC Scout	12.00	2.40	MSL-5389NIMH
	PC1 90303	12.00	2.50	GL8108
	Qube	0.00	0.00	OEM ONLY
	TEK 400	4.80	7.00	M400
	Tek 400	4.80	7.00	M400
	TEK 414 413A	4.80	7.00	MS413A
	UCW Monitor	8.40	0.10	9V-NIMH
	UCW Monitor	3.00	1.30	DL123AB
	VITATEK 500	16.00	2.50	GL1625
Spiro Medix				
	3300	1.20	0.60	MSL-S3300
Sscor				
	2314	12.00	2.80	MSL-5052
Stanley				
	Card Reader Alkaline	6.00		MSL-6219
Stereo-Optical				
	Retnascope 360-H11	3.60	0.30	MSL-5738
Steris				
	3080 Control	12.00	2.80	L12V2.4F1-06
	3080 Motor	12.00	28.00	L12V28B1-15

COMPANY	MODEL	VOLTS	AH	BATTERY
Steris				
	3085 (B00815-1)	12.00	28.00	L12V28B1-15
	3085 (B10998)	12.00	2.80	L12V2.4F1-06
	4085	12.00	28.00	NPX-100RFR
	C-Max	12.00	26.00	L12V26B1-16
	ESC Chair	12.00	1.20	L12V1.2F1-04
Stochard				
	S5	12.00	18.00	L12V18B1-14
Stryker				
	3000	3.60	0.50	MSL-5774
	5920-010-091	0.00	0.00	OEM ONLY
	Drill	0.00	0.00	OEM ONLY
	Smartpump Dual Channel	0.00	0.00	OEM ONLY
	System 6 Drill	0.00	0.00	OEM ONLY
	T4 400-650	6.00	0.38	MSL-400-650
Sunrise Medical				
	7304A Suction	12.00	1.30	L12V1.3F1-04
	7305P-D	12.00	2.50	MSL-5975
Suntech				
	247	6.00	5.00	L06V4.5F1-01
Surgidyne				
	247B	6.00	5.00	L06V4.5F1-01
	350-21	7.20	1.40	MSL-5243
	350-2I	7.20	1.40	MSL-5243
	500-22E	12.00	7.20	L12V7.2F1-04
Sylvan				
	All Models	3.60	0.15	MEMPAK3.6V
	Pediascan	0.00	0.00	OEM ONLY
Symbol				
	82-67705-01	3.60	0.70	BCS-LS4278
TDI Circuit Guard				
	CGV	3.60	0.15	MEMPAK3.6V
Technidyne				
	Hemochron401	7.20	1.30	MSL-5479
Thermo Scientific				
	REB2304	12.00	1.30	L12V1.3F1-04
	TSX505SA	3.00	1.40	CR2/3A
Tollos				
	Medi-Lifter III	12.00	33.00	L12V33B1-15
	Titan-X (Pat Lift)	12.00	33.00	L12V33B1-15
Tono-Pen				
	XL	3.00		CR1/3N
TrippLite				
	BP240V10RT3U*	12.00	9.00	UP-VW1245P1
	Smart500RT1U	12.00	7.00	MSL-2375/AP
	Smart500RT1U (Req 2)	6.00	7.20	L06V7.2F1-02
Tyco				
	Kangaroo Control	7.20	1.60	MSL-5070
Varidyne				

COMPANY	MODEL	VOLTS	AH	BATTERY
Varidyne	V350 *	12.00	7.20	L12V7.2F1-04
Verathon	0400-0100	0.00	0.00	OEM ONLY
	0400-0127	0.00	0.00	OEM ONLY
	0570-0338	0.00	0.00	OEM ONLY
	0800-0421	0.00	0.00	OEM ONLY
	BVI3000	0.00	0.00	OEM ONLY
	BVI5000	0.00	0.00	OEM ONLY
	Glidescope 0400-0064	0.00	0.00	OEM ONLY
	Glidescope 0400-0102	0.00	0.00	OEM ONLY
	Glidescope CZ192LIP	0.00	0.00	OEM ONLY
	Glidescope GVL	0.00	0.00	OEM ONLY
	KNBNK513475	0.00	0.00	OEM ONLY
Viasys	Avea	24.00	4.50	MSL-6068
	Avea External	12.00	12.00	L12V12F2-04
Vickers	4000	4.00	1.00	L04V1.0C1-17
Vital Air	LTV1000	12.00	12.00	L12V12F2-04
	LTV1000	6.00	5.00	L06V4.5F1-01
Waverly Glenn	C-450*	12.00	4.50	L12V5.0F1-03
	C625*	12.00	4.50	L12V5.0F1-03
	F600	24.00	2.90	MSL-PL-WB-00200
Welch Allyn	008-0125-00	8.00	5.40	MSL-5503
	010-0230-03	0.00	0.00	OEM ONLY
	105631	0.00	0.00	OEM ONLY
	105632	0.00	0.00	OEM ONLY
	300 Series	6.00	4.50	MSL-6065
	407560	6.00	4.50	MSL-6065
	4200	6.00	4.50	MSL-5077
	4200-84	6.00	4.50	MSL-6039
	42NOB	6.00	4.50	MSL-6039
	42NTB Li-Ion	0.00	0.00	OEM ONLY
	4500-84	6.00	6.00	MSL-6131
	45NTO	6.00	6.00	MSL-6131
	501-0015-01	6.00	4.50	MSL-6065
	52000	6.00	4.50	MSL-5077
	53000	6.00	4.50	MSL-6065
	53NTO	6.00	4.50	MSL-6065
	53NTP	6.00	4.50	MSL-6065
	6000	0.00	0.00	OEM ONLY
	61001	3.60	0.10	MS72500
	62000	6.00	7.20	L06V7.2F1-02
	71150	3.50	0.70	MS722
	71670	3.50	0.70	MS722

COMPANY	MODEL	VOLTS	AH	BATTERY
Welch Allyn	72000	2.50	1.00	MS720
	72100	2.50	0.70	MS721
	72200	0.00	0.00	OEM ONLY
	72240	4.80	1.80	MS72240
	72250	3.60	1.80	MSL-6045
	72270			MS72270
	72300	3.60	0.70	MS72300
	72500	3.60	0.10	MS72500
	72600	2.40	0.60	MSL-5182
	72700	7.20	0.10	MSL-5511
	75200	3.60	1.80	MSL-6045
	Atlas	6.00	7.20	L06V7.2F1-02
	Atlas 6200	6.00	7.20	MSL-5053
	BATT11	0.00	0.00	OEM ONLY
	BATT22	0.00	0.00	OEM ONLY
	BATT69	0.00	0.00	OEM ONLY
	BATT99	10.80	6.75	MSL-BATT99
	Connex	10.80	6.75	MSL-BATT99
	CP10	12.00	2.00	L12V2.0F1-09
	CP100	6.00	4.50	MSL-5077
	CP200	6.00	4.50	MSL-5077
	Life Sign	6.00	5.00	L06V4.5F1-01
	Lumi View, 72250	3.60	1.80	MSL-6045
	ProBP 2400	0.00	0.00	OEM ONLY
	ProBP3400	0.00	0.00	OEM ONLY
	Propaq 242*	8.00	2.70	L08V2.5F1-04
	Propaq CS	8.00	5.40	MSL-5503
	Propaq LT	0.00	0.00	OEM ONLY
	Spot LXI	6.00	6.00	MSL-6131
	Spot420 VSM	6.00	4.50	MSL-6039
	Spot420 VSM	6.00	5.00	L06V4.5F1-01
	VSM6000 Series	10.80	6.75	MSL-BATT99
	VSN 42MTO	0.00	0.00	OEM ONLY
Wheel Co				
	Prevailer	12.00	30.00	L12V35B1-15
Zimmer				
	ATS1500	12.00	2.60	L12V2.6F1-02
	ATS2000	12.00	2.60	L12V2.6F1-02
	ATS3000	12.00	4.50	L12V5.0F1-03
	ATS4000	0.00	0.00	OEM ONLY
Zinc Air				
		1.40	0.00	EL123AP
Zoll				
	8000-0299-01	10.00	2.50	PD4420
	8019-0535-01	0.00	0.00	OEM ONLY
	AED Plus (req. 10)	3.00	1.30	DL123AB
	AED Pro	0.00	0.00	OEM ONLY
	E Series	0.00	0.00	OEM ONLY

COMPANY	MODEL	VOLTS	AH	BATTERY
Zoll	M Series	10.00	2.50	PD4420
	NTP1000	12.00	4.50	MS461
	PD1200	12.00	2.50	90-M1200PD
	PD1400	10.00	2.50	PD4420
	R Series	0.00	0.00	OEM ONLY