

# Specialty Bulbs

*Medset Specialties Ltd*  
*3334 Mainway Drive*  
*Burlington, On L7M 1A7*

**Tel: 800-331-3187**

**Fax: 905-335-4456**

Health Canada Establishment License # 1810

**Email: [sales@medset.ca](mailto:sales@medset.ca)**

**Web: <http://medset.ca>**

COMPANY	PART#/MODEL	VOLTS	WATTS	BULB
3M				
	00-1912	0.0	0	OEM ONLY
	2000AK	12.0	50	EXZ
	788045 45715	12.0	75	JCRM12V75W/H
ACMI/Circon				
	FCB96	21.0	150	EJA
	Gyrus MV9090			CL1583
	HEADLIGHT	21.0	150	EJA
	MV9082L	45.0	250	HTI250W/32
	MV9086L	0.0	300	90-9300BF
	XLS-3000			BT-PE300C10F
Aimes				
	CLINATEK 200	11.0	30	EKZ
Air Shields				
	PT533	0.0	0	OEM ONLY
	PTM78-1,2,3,4	12.0	50	EXZ/CC
	PTS681	12.0	50	EXZ/PH
ALM Surgical				
	50067	23.0	100	ALM23V100W
	74000	24.0	120	74000
American Optical				
	PROJECTO CHRT	125.0	75	CAX BT
Amersham Phar				
	Ultraspec 1000	12.0	20	64425
Amsco				
	129359-934	22.0	200	P129359-934
	764-319-682	33.0	235	P129249-001
	AS4810	22.0	200	P129359-934
	EXAM 10	120.0	150	ESP
	GEMINI AS4705	0.0	0	OEM ONLY
	JD1073	120.0	150	ESP
	OR LIGHTS	10.0	150	ESP-WK
	Q235T4/3	33.0	235	P129249-001
	QUANTUM	22.0	220	P129362-228
	SOCKET	0.0	0	SOC-P129249-00
AO				
	1810	7.0	17	1493
Bausch & Lomb				
	31-13-21	6.0	20	ESB
	71-71-84	120.0	15	71-71-84
	71-71-84	120.0	15	XXGE15S111412

COMPANY	PART#/MODEL	VOLTS	WATTS	BULB
Bausch & Lomb	LEMSOMETER	120.0	15	71-71-84
Baxter	7100	0.0	0	09500
Beckman	596791	0.0	0	SCD753PF
	MODEL 25	0.0	0	25-001H
Bell & Howell	PROJECTOR	24.0	250	ELC
Berchtold/Marti	C570	24.0	150	BV1413
	C571D	24.0	150	BV1413
	C950C	12.0	50	ATM36
	Chromophare D300	12.0	50	64602
	CZ7955	12.0	50	ATM36
	CZ904-24	22.8	50	018566US
	CZ908-22	22.8	150	CZ908-22
	CZ908-24	24.0	150	BV1413
	CZ908-5004		70	908-5004
	CZ-909-22	24.0	150	64465U
	CZ940-12	12.0	35	JC12V35W
	CZ955-12	12.0	50	64602
	JC12V35W	12.0	35	JC12V35W
BFW Inc	MAXENON	0.0	300	90-9300BF
BLV	ENDOSCOPE	24.0	100	BLV115610
Boehm	SL40	4.0	0	03100
Burton	1006926	120.0	250	MODIFIED UST
	9110 (0315010)	115.0	0	FHTSBLB
	CoolSpot 13410	20.0	150	DDL-US
	LAMP	120.0	100	100G161/229SC
	SPOTLIGHT	120.0	250	ENH
Cabot	2974	45.0	250	HTI250W/32DYO
Canon	BFO-2241-020-0	0.0	0	BFO 2241 020-0
	BH-33015-000	0.0	0	BH-33015-000
Castle		120.0	150	Q150CL/DC
		120.0	150	ETC
	110156	24.0	100	BLV115610
	150P25/10	120.0	150	150P25/10
	150P25/2SB	120.0	150	P25 CLSB
	2400	0.0	0	ESP-SOC
	2400 SERIES	120.0	150	ESP
	JD1073	130.0	150	EYJD130V150WB

COMPANY	PART#/MODEL	VOLTS	WATTS	BULB
Castle	JD1073	120.0	150	ESP
	K448	24.0	60	K448
	LS40	120.0	75	64484/FR
	MDT 516767	24.0	60	K448
	OR LIGHT	120.0	150	ETC
	Q150CL/CL	120.0	150	ETC
Chick Light	74000	24.0	120	74000
Coherent	SLIP LAMP	12.0	30	64260
Coltene-Wholed	COLTOLUX 75	12.0	75	64617
Concept		45.0	250	HTI250W/32CO
	8426	45.0	250	HTI250W/32DYO
	LIS8430	0.0	300	90-9300BF
Contec	LS6000	0.0	0	XBO-R100W/45-
Copeland-Optec	OPTEC 360	2.0	0	BT360125
Cuda	M2-150	21.0	150	DNF
	MX2-300	45.0	250	HTI250W/32CUD
	XLS180	0.0	0	XBOR180W-45C
	XLS300		300	90-8600PE
Curtis	TIMER	0.0	0	TM3
Dentsply	644751	10.0	52	13298
Dialux	20EB	6.0	20	ESB
Downs Surgical	QE-020-01	12.0	20	GBE
Dr. Mach		22.8	50	018566US
	130	22.8	50	018566/H-US
	67100201	22.8	50	018566/H-US
Draeger	AnestAnesthesia cart	12.0	20	FTD
	DUL18W31-830	61.0	18	DUL18W31-830
	H600W	12.0	35	FMW
	MU03690	12.0	50	EXN
	PT4000 (White)	61.0	18	FT18DL/830/EC
Dyonics	2974	45.0	250	HTI250W/32DYO
	300XL	0.0	300	90-9300BF
	310-752	45.0	250	HTI250W/32DYO
Emerg Light				

COMPANY	PART#/MODEL	VOLTS	WATTS	BULB
Emerg Light		12.0	0	908
Fiber Optics ME	INFINITY	6.0	10	ESA/FHD
	TRANSILLUMINAT	6.0	10	64410
Fiche		14.0	25	FLT
Florescent		0.0	0	F20T12-D
Fujinon	EXP-2200		150	PE150AF
	FLX-2000	0.0	300	90-9300BF
	FLX-300A	0.0	300	90-9300BF
	LX-300-F2	0.0	300	90-9300BF
GE/Marquette		82.0	410	FXL
	Giraffe/Panda	12.0	35	M199/FTH
	M54FST	12.0	20	FST-US
General Marc	GEMINI 300	35.0	300	EZG
	GEMINI 300-16	38.0	300	EZM
General Radio	1531-P1	6.0	6	157
Geraeus	Overhead Lamp	22.8	50	018566/H-US
Gilras	GR-7 Slit Lamp	6.0	20	64251HLX
Gilway	TECHLOOP	12.0	20	64425
Global	74000	24.0	120	74000
	ECA002	24.0	120	74000
Goldman	PERIMETER	6.0	5	5006
	PERIMETER	6.0	27	940-750
Haag Streit		6.0	26	BT366
	900	6.0	27	BT3050F
	900 930	6.0	26	BT366
	BD900	12.0	30	HS1200559
	BQ900	6.0	27	BT3050F
	HS0900.1551	6.0	1	MC1009P
Hanalux	018366	22.8	75	018366-BT
	018366	22.8	75	018366US
	018566	23.0	50	018566-DF
	018566	22.8	50	018566-HX
	018769	23.0	40	DF94-60-31
	2000	22.8	50	018566/H-US

COMPANY	PART#/MODEL	VOLTS	WATTS	BULB
Hanalux				
	560 53 010	24.0	50	CL70549
	561-15282	22.8	50	018566/H-US
	64650	23.0	50	64650
	Blue 30	24.0	50	CL70549
	Blue 90	22.8	90	53198
Haskel				
		14.0	2	MIN-161
Healthco				
	ARISTOCRAT OTP	12.0	55	64151
Heine				
		3.0	0	X-01-88-059
		2.5		X-01-88-105
		4.0	0	X-02-88-078
		4.0	0	X-02-88-086
		6.0	0	X-04-88-060
		6.0	0	X-04-88-068
	EN10-S	4.0	0	X-02-88-046
	HLX 046	4.0	0	HEX 02-88-046
	OPT/SCOPE	4.0	0	X-02-88-070
	OTOSCOPE	3.5	0	X-02-88-044
	X-01.88.035	2.5	0	06000
	X-01.88.056	2.5		X-01-88-056
	XHL060	6.0	0	X-04-88-060
Heraus				
	2000	22.8	50	018566/H-US
Hewlett-Packard				
		14.0	0	7382
	3964	0.0	0	PRE-ALIGNED
Hiber-Teck				
	SHIMADZU	12.0	50	64610
InFocus				
	X2/X3	0.0	0	OEM ONLY
Isolux				
	EX300C	0.0	300	90-9300BF
Jena				
	681-34/3	6.0	15	NA67271
Keeler				
	1012-P-7003	6.0	10	1012-P-7003-CL
	VANTAGE	6.0	10	1012-P-7003
	VANTAGE (alternate)	6.0	10	1012-P-7003-CL
Kodak				
	EKTAGRAPHIC 3	82.0	300	EXR
	PROJECTOR	82.0	300	FHS
Kowa				
	SL-5	12.0	20	BT12V20W
Laborie/Lumatec				
		0.0	300	90-9300BF
		21.0	150	EJA

COMPANY	PART#/MODEL	VOLTS	WATTS	BULB
Laerdal				
	SUCTION	5.0	0	7333
Leica				
	1038 4643	12.0	50	HLX4643
	115003	12.0	30	64261
	3040-004	0.0	50	HBO50W/L2
	3040-0090	0.0	50	HBO50W/L1
	3040-0150	6.0	15	8018
	3040-0280	6.0	20	FHE/ESB
	3040-0290	6.0	10	5001
	8031844	6.0	15	8018
	ACT2000	6.0	9	13529
Leisengang				
	1DSUL	12.0	20	64425
Leitz				
	1008G	6.0	27	BT3050F
	70302	12.0	60	348
	HS360	6.0	27	BT3050F
	Laborlux K	6.0	20	ESB
	LABORLUX S	12.0	100	FCR
	MICROSCOPE	6.0	10	ESA/FHD
Linvatek				
				BT-PE300C10F
	8425/8426	45.0	250	HTI250W/32
	AT900	6.0	20	ESB
	C3141	0.0	0	XBOR180W-45C
	LIS8430	0.0	300	90-9300BF
	LS7500	0.0	300	XBOR300W/60C
	LS7700			57980
Luxo				
	LMP-1	120.0	100	CEA
Luxtec				
	1300A	21.0	150	DNF
	9300XSP	0.0	300	90-9300BF
	Series ACO	15.0	150	EFR
	Y1144	0.0	0	91-150F
MacBeth				
	SPECTRA	120.0	750	EJG
Mallincrodt				
	12866-5606	12.0	10	6411
MARCO				
	41340	6.0	26	BT366
	CP690E	6.0	20	ESB
Martin				
	908-24	24.0	150	BV1413
	ATM 36	24.0	250	CZ907-24
	B18CZ-906-24	24.0	150	BV1413
	C300	22.8	50	018566US
	C570	24.0	150	BV1413

COMPANY	PART#/MODEL	VOLTS	WATTS	BULB
Mavig	Uniflex	23.0	50	64650
MDT/Castle	SERIES 2000	24.0	60	K448
Medical Illimin	0001256	12.0	90	64458
	011617	12.0	50	64440AX-OS
	11617	24.0	150	BV1413
	Centurion 021513	12.0	90	64458
Medicana	561-15282	22.8	50	018566/H-US
	HEUREUS	22.8	50	018566/H-US
Metler	P1200	6.0	10	657
Midmark	6600	24.0	250	64654
	77001	0.0	0	OEM ONLY
Miles	CLINATEC 200+	0.0	0	EKZ-CLIN200
Narva	67271	6.0	15	NA67271
NEC	440	0.0	160	UMRET400LT
Nicatron	044	3.5	0	X-02-88-044
	044	4.0	0	BOX0288044
Nicolet	VOYAGEUR	0.0	42	FGQG33-1.75
Nikon	77458	115.0	30	77458
	77903	6.0	15	77903
	77904	6.0	5	77904
	Eclipse E400	6.0	30	5761
	Opti-Phot	12.0	50	64610
	R400	6.0	0	328
Ohmeda	3100 Infant Warmer	12.0	50	EXN-PH
	3300 Exam Light	12.0	50	EXN-PH
	6600-0160-800	21.0	150	EKE-USHIO
	6600-0235-850	0.0	0	OEM ONLY
	6600-0262-200	20.0	150	DDL-US
	6600-0680-200	20.0	150	6600-0680-200
	Giraffe	0.0	0	OEM ONLY
	Panda Warmer Alarm	120.0	6	0690-2100-315
	PHOTOTHERAPY	120.0	20	F20T12BB
	PHOTOTHERAPY	0.0	0	OEM ONLY
Olympic	51418	0.0	0	OEM ONLY
	Bili-Lite-Pad 91	21.0	150	EKE-USHIO

COMPANY	PART#/MODEL	VOLTS	WATTS	BULB
Olympus		6.0	10	8B104
		0.0	300	90-9300BF
	217073	0.0	0	TM3
	23287	15.0	150	EFR
	3-6V2AGB-4	0.0	0	OEM ONLY
	58107	0.0	300	90-9300BF
	58111	0.0	300	90-9300BF
	7850	0.0	0	OEM ONLY
	8C101	6.0	5	77904
	8C401	6.0	30	5761
	8-C410	12.0	50	EXN
	A1/232	15.0	150	EFR
	BH2	0.0	100	HBO100W/2-US
	BH2	0.0	100	HBO100W/2
	BX41	6.0	30	5761
	CL1584	0.0	0	OEM ONLY
	CLE-10	15.0	150	EFR
	CLE3	15.0	150	EFR
	CLE-44	15.0	150	EFR
	CLE-ALL	15.0	150	EFR
	CLH-SC	15.0	150	EFR
	CLK-3	15.0	150	EFR
	CLK4	15.0	150	EFR
	CLV 140	0.0	300	90-9300BF
	CLV-10	0.0	300	90-9300BF
	CLV-160	0.0	300	90-9300BF
	CLV-180	0.0	300	90-9300BF
	CLV-190		300	90-9300ES
	CLV-290SL		300	90-9300ES
	CLV-S190	0.0	300	90-9300BF
	CLV-S40	0.0	300	90-9300BF
	CLV-U20	0.0	300	90-9300BF
	CLV-U40	0.0	300	90-9300BF
	CV1	0.0	300	90-9300BF
	GB4	0.0	0	OEM ONLY
	GB-4BH2	0.0	0	OEM ONLY
	I75LS1	0.0	300	90-9300BF
	K10001866	15.0	150	EFR-OS
	MAJ 525			MAJ525
	MAJ Series	0.0	0	OEM ONLY
	MAJ1817		300	90-9300ES
	MB-361	0.0	300	90-9300BF
	MD-151	15.0	150	EFR-OS
	MD631	0.0	300	90-9300BF
	MICROSCOPE	8.0	40	8B101
	MICROSCOPE	6.0	5	8G101
	MICROSCOPE	8.0	30	8C102
	OCS 500	15.0	150	EFR-OS



COMPANY	PART#/MODEL	VOLTS	WATTS	BULB
Olympus	POINTER	6.0	10	8G102
	SPOTLIGHT	12.0	50	EXT
	SZ61TR	12.0	75	JCRM12V75W/H
	TIMER	0.0	0	TM2
	Y1052	0.0	300	90-9300BF
	Y1064	0.0	300	90-9300BF
	Y1064SP	0.0	300	90-9300BF
	Y1176	0.0	300	90-9300BF
Ortholux				
		6.0	5	8100
Osram	RING DAYLIGHT	120.0	32	FC12T9/D/RS
	XPR120W/45CVI			XPR120W/45CVI
Ozinator				
		0.0	0	GPH287T5PH
Penlon				
		6.0	0	55094
	04700	3.0	0	55092
	LARYNGOSCOPE	3.0	0	55091
Pentax				
		0.0	500	LX500A
	5000i	0.0	0	OEM ONLY
	7000i	0.0	0	OEM ONLY
	ATM3300	0.0	300	90-9300BF
	BS-H2			N/A
	EMP3000	0.0	300	90-9300BF
	EPK1000	13.0	100	OLX20
	EPK700	0.0	0	OEM ONLY
	EPK-I		300	CL1576
	LH-150	21.0	150	DNF
	LNS-BP310			N/A
	LX-300	0.0	300	90-9300BF
	LX-75F	0.0	75	XBO75W/2-OFR
	LX75FSY518	0.0	300	90-9300BF
	OL-X4	0.0	300	90-9300BF
	OLX-4	0.0	300	90-9300BF
Perkin Elmer				
	C0060501	9.0	35	BWD
Philips				
		6.0	15	8018
		6.0	5	5006
		6.0	20	ESB
		15.0	150	EFR
	6423	15.0	150	EFR
	bSureLC3141	0.0	0	LCA-3118
	MICROSCOPE	6.0	15	13347C
	PL-L36W/01			N/A
Radiag				

COMPANY	PART#/MODEL	VOLTS	WATTS	BULB
Radiag		11.0	30	DZA
Radx		120.0	0	RADX
	ODAMEX	120.0	16	FLU16/D
Reichert	PROJECTOR	6.0	20	ESB
Respironics	ILS-100 / ILS-101	6.0	10	64410
	Infinity	6.0	10	64410
Rusch	GreenSpec Laryng.	2.5	2	5630
Rush-Pilling	52119	15.0	150	EFR
	521285	15.0	150	EFR
	529180	0.0	0	XBOR180W-45C
	LUMINATOR	21.0	150	DNF
	LUMINATOR 2X	15.0	150	EFR
S&S Xray	118	120.0	75	PH/211
Sartorius	69/0315	6.0	6	818-47-43
Shemadzu	UV160	0.0	0	DS-350U
Shin-Nippon	SL-40	6.0	20	ESB
	SL-40	6.0	1	1768
Siemens		24.0	150	64642-OS
	56115282	22.8	50	018566/H-US
	MAMOGRAPHY	12.0	50	64440AX
	MAMOGRAPHY	12.0	50	64440
Sigma	ANALYZER	14.0	85	13288
Simms		24.0	50	64445
Skylux	ST-HCR153S	24.0	40	XXSTHCR1503S
Skytron		24.0	40	A1-010-28-1
	294-B1-010-28	24.0	50	BTH24501
	90 & 58	24.0	40	HCR15035
	98	0.0	0	OEM ONLY
	A1-010-18	0.0	0	OEM ONLY
	A1-010-28-1	24.0	40	A1-010-28-1
	B1-010-28	24.0	50	BTH24501
	B1-011-32	24.0	100	DKK
	B4-010-29	24.0	75	B4-010-29
	B5-010-04	0.0	0	B5-010-04

COMPANY	PART#/MODEL	VOLTS	WATTS	BULB
Skytron	EX62	24.0	50	BTH24501
	HA60AR	24.0	100	DKK
	HALOGEN	24.0	100	DKK
	HOR-150-3S	0.0	0	OEM ONLY
	SH42	24.0	40	A1-010-28-1
	SH62	24.0	60	A1-420-010
	STELLAR	24.0	50	BTH24501
	WAND	0.0	0	B5-010-04
Socket				
	UNIVERSAL	0.0	0	32400
Socket EJM				
		0.0	0	QLV-1 (FRAME)
Spectra Electri				
	MICROSCOPE	30.0	80	ELB
Spectronics				
	BLE480B	0.0	0	F6T5BLB
Steri-Lite				
		0.0	0	GT10T51/2L
Steris				
	Harmony 500/700	20.0	180	P093926-047
	P084079-033	120.0	150	ESP
Stortz				
	201315-20	0.0	175	90-9175BF
	201320-20	0.0	175	90-9175BF
	201321-20	0.0	175	90-9175BF
	201330-28	0.0	300	90-9300BF
	201331-20	0.0	300	90-9300BF
	201340-01	0.0	300	90-9300BF
	264H	0.0	300	90-9300BF
	482CA	15.0	150	BRJ
	487XE	0.0	300	90-9300BF
	610C	0.0	300	90-9300BF
	611	0.0	300	90-9300BF
	611C	0.0	300	90-9300BF
	615C	0.0	175	90-9175BF
	615XE	0.0	300	90-9300BF
	615XES	0.0	300	90-9300BF
	630	0.0	300	90-9300BF
	7748	24.0	250	64655-OS
	LX175-F	0.0	175	90-9175BF
	MILLENNIUM	0.0	0	WAM21E00S-00
	N009	6.0	6	157
	N1002	6.0	6	157
	Nova 300	0.0	300	90-9300BF
	Nova175	0.0	175	90-9175BF
	ST103	15.0	150	BRJ
	Y1441	0.0	175	90-9175BF
Stryker				

COMPANY	PART#/MODEL	VOLTS	WATTS	BULB
Stryker		45.0	250	HTI250W/32STR
		45.0	250	HTI250W/22STR
	2112	28.0	70	JCD28V70W
	2188-1503-58	28.0	70	JCD28V70W
	220-180-300	45.0	250	HTI250W/22S-ST
	220-191-000			BT-PE300C10F
	5000	45.0	250	HTI250W/22S-ST
	BRACKET	0.0	0	MSL-Q5000B
	QUANTUM 3000	45.0	250	HTI250W/32STR
	QUANTUM 4000	45.0	250	HTI250W/22STR
	QUANTUM 5000	45.0	250	HTI250W/22S-ST
	QUANTUM 5000	45.0	250	HTI250W/22
	QUANTUM 6000	0.0	300	90-9300BF
	X6000	0.0	300	90-9300BF
	X7000			BT-PE300C10F
	X8000			ASK FOR DETAILS
	X8000 (No Timer)			HY-CL1567
	X8000 Bulb			HY-CL1495
Sylvania				
	CIRCLINE	120.0	0	F15W/154-T8
	PROJECTOR	19.0	80	DDM
Topcon				
	SL-2D	6.0	20	403-402-0700
	SL-8Z	12.0	30	44680-25700
	SL-D7	12.0	30	44680-25700
TRAFOS				
		6.0	0	X-04-88-060
	E7	6.0	0	HEX0488060
Tyco/Simplex				
		12.0	6	CFM12V6.2M
Vanguard				
		120.0	500	BCK
Volpi				
	4000-1	15.0	90	EPV
Voltarc				
		120.0	0	FS72T12/UVB/H
Weck				
	173301	0.0	300	90-9300BF
	521118	21.0	150	DNF
	521311	0.0	0	XBO-R100W/45
	801A202	6.0	30	WKEI738
	90-6306	0.0	300	90-9300BF
	FRX-2000	0.0	300	90-9300BF
	WECKLITE300	0.0	300	90-9300BF
Wehmer				
	HEADLAMP	6.0	10	ENT
Welch Allyn				
	00200	3.0	0	00200

COMPANY	PART#/MODEL	VOLTS	WATTS	BULB
Welch Allyn				
	00900	3.0	0	00900
	01100	3.0	0	04700
	02000	0.0	0	02000
	02500	6.0	6	157
	02600	0.0	0	02600
	02700	3.0	0	03400
	03000	4.0	0	03000
	03000	4.0	0	03000-WA
	03100	4.0	0	03100-WA
	03100	4.0	0	03100-Z
	03100	4.0	0	03100
	03300	3.0	0	03300
	03400	3.0	0	03400
	03800	3.0	0	03400
	03900	0.0	0	03900
	04000	6.0	0	07800
	04100	0.0	0	04100
	04200	12.0	35	EPN
	04400	3.0	0	04400
	04700	3.0	0	04700
	04800	3.0	0	04800
	04900	4.0	0	04900
	06000	2.5	0	06000
	06300	12.0	20	64425AX
	06400	0.0	0	ESX/FG
	07600	0.0	0	ESX/FG
	07800	6.0	0	07800
	08100	12.0	20	64425
	08200	0.0	0	08200
	08500	0.0	0	08500
	08800	0.0	0	08800
	09500	0.0	0	09500
	20T3Q/CL	12.0	20	64425
	25020	4.0	0	03100
	48100	12.0	20	64425
	48200 Lite Pipe	12.0	35	EPN
	69601	2.5	0	06000
	74700OP	4.0	0	03000
	78810	0.0	0	08800
	901079	4.0	0	03100
	HPX060	3.0	0	04800
	LS135	12.0	35	EPN
	LS150	12.0	50	04450
	LS200	0.0	0	ESX/FG
	POCKET SC OPTL	3.0	0	03400
	POCKET SC OPTO	3.0	0	03400
Wilde				
	1038-4643	12.0	50	HLX4643

COMPANY	PART#/MODEL	VOLTS	WATTS	BULB
Wilde	13846430	12.0	50	HLX4643
	384643	0.0	0	OEM ONLY
	M-20, M-12	6.0	20	70249
	M655	12.0	50	HLX4643
	M655	12.0	50	Wild50
Wolf	4200	24.0	150	64640-OS
	4251	24.0	250	ELC
	D4002U	21.0	150	EJM
	D4008U	21.0	150	EJM
	LAPROSCOPE	55.0	400	HTI400/24
XRite	Densitometer301			F6T5-D
Zeiss	234-755	12.0	50	64609
	301	0.0	300	90-9300BF
	301387	6.0	20	64251HLX
	304953-9030	0.0	300	90-9300BF
	310752	45.0	250	HTI250W/32DYO
	380018-1740	6.0	15	77Z
	380029-9020	6.0	15	8017
	380059-1660	12.0	100	FCR
	380079-9030	12.0	75	EFN
	380079-9040	12.0	100	EFP
	380079-9230	12.0	100	380079-9230
	380079-9240	22.8	50	018566
	380079-9550	12.0	50	64440AX
	38-01-20	12.0	30	64260
	380-120-7030	6.0	10	64222
	380120-7040	12.0	30	64260
	380177	6.0	15	13347W
	390153	6.0	25	BT53Z
	390158	6.0	30	BT58Z
	390186	6.0	50	BT86Z
	390413	110.0	50	390413
	407-309	0.0	0	XBOR180W-45C
	417053	15.0	150	EFR
	58102	12.0	100	EHE
	AXIO - A1	12.0	35	FRB
	AXIOVERT 25	6.0	25	HLWS-5A
	EYE PROJECTOR	12.0	100	EFP
	HXP120C			XPR120W/45CVI
	MICROSCOPE	120.0	6	S6
	MICROSCOPE	12.0	100	EFP
	MICROSCOPE	12.0	100	EHE
	OPMI	6.0	30	BT58Z
QG6AC	45.0	250	HTI250W/32DYO	
Slit Lamp	6.0	20	403-402-0700	

COMPANY

PART#/MODEL

VOLTS

WATTS

BULB

Zeiss

SUPERLUX

0.0

175

90-9175BF

SUPERLUX40

45.0

250

HTI250W/32

Y1153

0.0

300

90-9300BF